


## QPRO-μQPRO

### STANDARD CURVE PREPARATION



Calibration curve samples used to prepare BSA (Bovine Serum Albumine)  
Cod. PRD083,0010 BSA Standard

Vial	Volume of Diluent		Volume and Source of BSA		Final BSA Concentration	
	Qpro	μQpro	Qpro	μQpro	Qpro	μQpro
A	0 μL	4.5 mL	300 μL of Stock 2.0 μg/mL	0.5 mL of Stock 200 μg/mL	2.0 μg/mL	200 μg/mL
B	125 μL	8.0 mL	375 μL of Stock 2.0 μg/mL	2.0 mL of vial A dilution	1.5 μg/mL	40 μg/mL
C	325 μL	4.0 mL	325 μL of Stock 2.0 μg/mL	4.0 mL of vial B dilution	1.0 μg/mL	20 μg/mL
D	175 μL	4.0 mL	175 μL of vial B dilution	4.0 mL of vial C dilution	750 μg/mL	10 μg/mL
E	325 μL	4.0 mL	325 μL of vial C dilution	4.0 mL of vial D dilution	500 μg/mL	5 μg/mL
F	325 μL	4.0 mL	325 μL of vial E dilution	4.0 mL of vial E dilution	250 μg/mL	2.5 μg/mL
G	325 μL	4.8 mL	325 μL of vial F dilution	3.2 mL of vial F dilution	125 μg/mL	1 μg/mL
H	400 μL	4.0 mL	100 μL of vial G dilution	4.0 mL of vial G dilution	25 μg/mL	0.5 μg/mL
I	400 μL	8.0 mL	0 μL	0 mL	0 μg/mL = Blank	0 μg/mL = Blank

### Related products



Downloads: <http://www.cyanagen.com/downloads/product-manuals#family-1>

**CYANAGEN**

Reagents for Molecular Biology

Cyanagen s.r.l. T +39 051 534063  
Via Stradelli Guelfi, 40/C info@cyanagen.com  
40138 Bologna (ITALY) www.cyanagen.com



**CYANAGEN**

Reagents for Molecular Biology

THE RIGHT LIGHT

**QPRO-μQPRO**

protein quantitation  
Bicinchoninic Acid (BCA) kit

02

[www.cyanagen.com](http://www.cyanagen.com)

#### QPRO-μQPRO

##### TECHNICAL DESCRIPTION

**QPRO** and **μQPRO** are Bicinchoninic Acid (BCA) based protein assays for highly sensitive colorimetric analysis. The BCA method is fast and easy, providing enhanced flexibility.

The Bicinchoninic acid assays (BCA) rely on the production of  $\text{Cu}^+$  which is converted into a violet-coloured solution when reacting with Bicinchoninic acid.

**BSA** (Bovine Serum Albumine) is the protein used as main reference in biological tests (as Bradford, BCA, BCA assays, etc.).

##### STORAGE CONDITIONS

Store at RT.

#### QPRO-μQPRO

##### FEATURES

###### Sensitive

Broad linear working range with excellent sensitivity (20 to 2000  $\mu\text{g/mL}$  for QPRO and 0,5 to 20  $\mu\text{g/mL}$  for μQPRO).

###### Versatile

Can be used with tubes and Microplates.

###### Universal

Less protein-protein variation than dye-binding methods.

###### Robust

Unaffected by most ionic and nonionic detergents such as Triton™ X-100 and SDS (1%).

###### Easy-to-use

Much easier and faster than the classical Lowry method.

#### QPRO-μQPRO

##### QUICK START PROTOCOL

Pipette standard and unknown sample into appropriately labelled tubes

- 0,1mL for QPRO
- 1 mL for μQPRO.

Add the Working Solution and mix

- 2 mL (Reagent A + Reagent B, 50:1 vol. ratio) for QPRO
- 1 mL (Reagent A + Reagent B + Reagent C, 25:25:1 vol. ratio) for μQPRO.

Cover tubes and incubate

- at 37°C for 30min or at RT for 2h for QPRO
- at 60°C for 60min for μQPRO.

Cool down to RT and measure absorbance at 562nm using a spectrophotometer and use water as a blank.

Use the standard curve to determine the protein concentration of the unknown sample.