


## RUBY GEL STAIN

### WARNING

The dye has optimal excitation at 302 and 470 nm, with an emission maximum at approximately 610nm.



### Related products



Downloads: <http://www.cyanagen.com/downloads/product-manuals#family-5>

**CYANAGEN**

Reagents for Molecular Biology

Cyanagen s.r.l. T +39 051 534063  
Via Stradelli Guelfi, 40/C info@cyanagen.com  
40138 Bologna (ITALY) www.cyanagen.com



**CYANAGEN**

Reagents for Molecular Biology

THE RIGHT LIGHT

# RUBY GEL STAIN

Fluorescent stain for  
protein detection in gels

03

[www.cyanagen.com](http://www.cyanagen.com)

## RUBY GEL STAIN

### TECHNICAL DESCRIPTION

**RUBY GEL STAIN** is a ready to use kit for rapid and sensitive protein staining in 1D and 2D SDS polyacrylamide electrophoresis gels.

It enables optimal visualization and quantitation of proteins with high contrast and the same sensitivity as silver staining without its drawbacks.

The staining procedure is a simple 220 minutes, three step protocol.

After staining, proteins can be removed from the gel and analyzed by mass spectrometry without interference from the stain.

### STORAGE CONDITIONS

Store at RT.

## RUBY GEL STAIN

### FEATURES

#### **Pure**

High purity of the dye: >98%.

#### **Optimized**

High signal to background ratio.

#### **Powerful**

Strong, uniform and reproducible signal from 0.2 ng to 10 ng protein.

#### **Time-saving**

Fast staining protocol (220 min).

#### **Complete**

Fixing and destaining solutions included in the kit.

## RUBY GEL STAIN

### QUICK START PROTOCOL

#### **Fixing Step (60 minutes)**

Incubate the gel in the Ruby Gel Stain Solution B for 1 hour using a rocker.

#### **Staining Step (60 minutes)**

Transfer the gel to a clean container and incubate it in the Ruby Gel Stain Solution A for 1 hour using a rocker. Wrap the staining container in aluminum foil to protect from direct light.

#### **Destaining Step (80 minutes)**

Transfer the gel to a clean container and destain it in Ruby Gel Stain Solution B for 80 minutes. Wrap the container in aluminum foil to protect from direct light.

#### **Washing (20 minutes)**

Transfer the gel to a clean container and equilibrate it for 20 minutes in deionized water and scan.