

CERTIFICATE OF CONFORMITY

IR-Blot 800 Goat anti-Mouse

Product Information

| | |
|--------------------------------|--|
| Product Number: | F1L108 |
| Lot Number: | AR2/17 |
| Description: | Secondary antibody optimized for multicolor infrared fluorescent Western blotting |
| Components: | Goat anti Mouse conjugate concentration 1mg/mL in PBS containing 0.01% sodium azide and 10mg/mL BSA after reconstitution |
| Storage: | +4°C |
| Stability/ Expiry Date: | 3 months after reconstitution at +4°C |
| Document Release Date: | 15/03/2018 |
| Document Revision: | 00 |

Conformance & Quality systems statement

This is to certify that this product conforms to Cyanagen specification.

The product was manufactured via a Quality System certified ISO9001 and tested in accordance with documented quality procedures and approved as a result of meeting the required specification.

Approved:



Quality Control Technician

Reviewed:



Production Manager

Cyanagen warrants that its products conform to the information contained in this and other Cyanagen publications. Purchaser must determine the suitability of the product for its particular use.