

CERTIFICATE OF CONFORMITY

Green Stain

Product Information


Product Number:	NAGS068
Lot Number:	IH29A-GR
Description:	Versatile and easy to use tool for the detection of nucleic acids in electrophoretic gels, based on post-electrophoresis gel staining, pre-electrophoresis gel staining and sample pre-staining.
Components:	10000 X in DMSO, 0.5 ml is sufficient for 50 mini-gels stain.
Storage:	-20°C, avoid repeated freeze thaw cycles.
Expiry Date:	05/2019
Document Release Date:	11/04/2018
Document Revision:	01

Conformance & Quality systems statement

This is to certify that this product conforms to Cyanagen specification.

The product was manufactured via a Quality System certified ISO9001 and tested in accordance with documented quality procedures and approved as a result of meeting the required specification.

Approved:



Quality Control Technician

Reviewed:



Production Manager

Cyanagen warrants that its products conform to the information contained in this and other Cyanagen publications. Purchaser must determine the suitability of the product for its particular use.