

CERTIFICATE OF CONFORMITY

QPRO-BCA kit Standard

Product Information


Product Number: PRTD1
Lot Number: MA09A-FB
Description: Bicinchoninic acid kit for protein determination
Components: Solution A: 1% BCA/tartrate in alkaline buffer, Lot MA09B-FB
Solution B: 4% CuSO₄ 5H₂O in water, Lot MA09C-FB
Storage: Room Temperature
Expiry Date: 01/2021
Document Release Date: 10/01/2019
Document Revision: 00

Conformance & Quality systems statement

This is to certify that this product conforms to Cyanagen specification.

The product was manufactured via a Quality System certified ISO9001 and tested in accordance with documented quality procedures and approved as a result of meeting the required specification.

Approved:



Quality Control Technician

Reviewed:



Production Manager

Cyanagen warrants that its products conform to the information contained in this and other Cyanagen publications. Purchaser must determine the suitability of the product for its particular use.

