

CERTIFICATE OF CONFORMITY

QPRO-BCA kit Standard

Product Information

Product Number:	PRTD1
Lot Number:	LH06A-LB
Description:	Bicinchoninic acid kit for protein determination
Components:	Solution A: 1% BCA/tartrate in alkaline buffer, Lot LH06B-LB Solution B: 4% CuSO4 5H2O in water, Lot LH06C-LB
Storage:	Room Temperature
Expiry Date:	08/2020
Document Release Date:	06/08/2018
Document Revision:	00

Conformance & Quality systems statement

This is to certify that this product conforms to Cyanagen specification.


The product was manufactured via a Quality System certified ISO9001 and tested in accordance with documented quality procedures and approved as a result of meeting the required specification.

Approved:



Quality Control Technician

Reviewed:



Production Manager

Cyanagen warrants that its products conform to the information contained in this and other Cyanagen publications. Purchaser must determine the suitability of the product for its particular use.