

CERTIFICATE OF CONFORMITY

SMART 488 – Antibody labeling kit

Product Information

Product Number:	F5A119N
Lot Number:	MB26A-XP
Description:	Kit for the fluorescent labelling (Abs 494 nm, Ex 520 nm) of antibody.
Components:	SMART 488 – Dye, Cod. F5A119N,DY, Lot MA29A-ET SMART – Elution buffer, Cod. SEB097,0250, Lot IF27A-NS SMART 488– Labeling buffer, Cod. SLB147,0003 Lot MB20A-ET SMART – Purification Resin, Cod. SPRO98, 0250, Lot IF27C-NS Empty plastic column
Storage:	SMART 488 – Dye at -20°C, other components at +4°C
Stability :	Components provide a kit performance stable up to 6 months from receipt, when properly stored.
Document Release Date:	26/02/2019
Document Revision:	00

Conformance & Quality systems statement

This is to certify that this product conforms to Cyanagen specification.


The product was manufactured via a Quality System certified ISO9001 and tested in accordance with documented quality procedures and approved as a result of meeting the required specification.

Approved:



Quality Control Technician

Reviewed:



Production Manager

Cyanagen warrants that its products conform to the information contained in this and other Cyanagen publications. Purchaser must determine the suitability of the product for its particular use.