

CERTIFICATE OF CONFORMITY

μQPRO-BCA kit Micro

Product Information

Product Number:	PRTD2
Lot Number:	IB06A-LB
Description:	Bicinchoninic acid kit for protein determination
Components:	Solution A: Tartrate in alkaline buffer, Lot IB06B-LB Solution B: 4% BCA in water, Lot IB06C-LB Solution C: 4% CuSO4 5H2O in water, Lot IB06D-LB
Storage:	Room Temperature
Expiry Date:	Components provide a kit performance stable up to 12 months from receipt, when properly stored
Document Release Date:	15/05/2019
Document Revision:	01


Conformance & Quality systems statement

This is to certify that this product conforms to Cyanagen specification.

The product was manufactured via a Quality System certified ISO9001 and tested in accordance with documented quality procedures and approved as a result of meeting the required specification.

Approved: 

Quality Control Technician

Reviewed: 

Production Manager

Cyanagen warrants that its products conform to the information contained in this and other Cyanagen publications. Purchaser must determine the suitability of the product for its particular use.