

**CYANAGEN**  
Reagents for Molecular Biology

**NEW**

THE RIGHT LIGHT

**BSA**

**Bovine Serum Albumine**

[www.cyanagen.com](http://www.cyanagen.com)

Cyanagen srl has a certified Quality System

ISO 9001:2008 QUALITY CERTIFIED

ISO 9001

BUREAU VERITAS  
Certification



Cyanagen s.r.l.  
Via Stradelli Guelfi, 40/C  
40138 Bologna (ITALY)

T +39 051 534063  
[info@cyanagen.com](mailto:info@cyanagen.com)  
[www.cyanagen.com](http://www.cyanagen.com)

# BSA

## TECHNICAL DESCRIPTION

**BSA** (Bovine Serum Albumine) is the protein used as main reference in biological tests (as Bradford, BCA, BCA assays, etc.).


Our BSA grants a high level of reliability, thanks to high purity and optimal characteristics of reproducibility in terms of concentration.

## FEATURES

- RELIABLE  
reliability in results
- STABLE  
extended shelf-life

Best results granted by using QPRO-BCA Kit Standar (PRTD1,0500) and  $\mu$ QPRO-BCA Kit Micro (PRTD2,0500).

BSA absorbance at 279 nm



product

**BSA Standard  
2,00 mg/mL**

cat #

PRTD083,0010

size

10 x 1 mL

**Cat #**