

## CERTIFICATE OF CONFORMITY

### μQPRO-BCA kit Micro

#### Product Information

Product Number:	PRTD2
Lot Number:	ND27A-CC
Description:	Bicinchoninic acid kit for protein determination
Components:	Solution A: Tartrate in alkaline buffer, Lot ND27B-CC Solution B: 4% BCA in water, Lot ND27C-CC Solution C: 4% CuSO <sub>4</sub> 5H <sub>2</sub> O in water, Lot ND27D-MZ
Storage:	Room Temperature
Stability/ Expiry Date:	Components provide a kit performance stable up to 12 months from receipt, when properly stored
Document Release Date:	27/04/2020
Document Revision:	00

#### Conformance & Quality systems statement

This document certifies that the product complies with the manufacturer specifications and its ISO 9001 Quality System.

The product was tested according to documented internal quality procedures and approved since it meets the required specification.

Approved:

  
Quality Control Technician

Reviewed:

  
Production Manager

Cyanagen warrants that its products conform to the information contained in this and other Cyanagen publications. Purchaser must determine the suitability of the product for its particular use.

**CYANAGEN SRL -**

C.F./P.IVA: 02360111203 iscritta nel Registro Imprese di Bologna n. REA: 433637

Capitale sociale sottoscritto e versato € 10.000,00.

"Società soggetta a direzione e coordinamento da parte di Finance & Technology Srl – art.2497 bis cc".