

CERTIFICATE OF CONFORMITY

Green Stain 10,000X in DMSO

Product Information


Product Code:	NAGS068
Lot Number:	OM08A-AR
Description:	Versatile and easy to use tool for the detection of nucleic acids in electrophoretic gels, based on post-electrophoresis gel staining, pre-electrophoresis gel staining and sample pre-staining
Components:	10000 X in DMSO, 0.5 ml is sufficient for 50 mini-gels stain.
Storage:	-20°C, avoid repeated freeze thaw cycles.
Stability/ Expiry Date:	07/2024
Document Release Date:	27/07/2023
Document Revision:	01

Conformance & Quality systems statement

This document certifies that the product complies with the manufacturer specifications and its ISO 9001 Quality System.

The product was tested according to documented internal quality procedures and approved since it meets the required specification.

Approved:



QC Manager

Reviewed:



Production Manager

Cyanagen warrants that its products conform to the information contained in this and other Cyanagen publications. Purchaser must determine the suitability of the product for its particular use.

MO 7.5.1.3 Rev04