

CERTIFICATE OF CONFORMITY

Q PRO – BCA Kit Standard

Product Information

Product Code:	PRTD1
Lot Number:	RC19A-AR
Description:	Bicinchoninic acid kit for protein determination
Components:	Solution A: 1% BCA/tartrate in alkaline buffer, Lot RC19B-AR Solution B: 4% CuSO ₄ 5H ₂ O in water, Lot RC19C-AR
Storage:	Room Temperature
Stability/ Expiry Date:	Components provide a kit performance stable up to 12 months from receipt, when properly stored
Document Release Date:	20/03/2024
Document Revision:	00

Conformance & Quality systems statement


This document certifies that the product complies with the manufacturer specifications and its ISO 9001 Quality System.

The product was tested according to documented internal quality procedures and approved since it meets the required specification.

Approved: 

QC Manager

Reviewed:



Production Manager

Cyanagen warrants that its products conform to the information contained in this and other Cyanagen publications. Purchaser must determine the suitability of the product for its particular use.

MO 7.5.1.3 Rev04