

CERTIFICATE OF CONFORMITY

STAR FLUOR 645 Antibody labeling kit

Product Information

Product Number:	F2H085N
Lot Number:	NM11D-XP
Description:	Kit for the fluorescent labelling (Abs 645 nm, Ex 665 nm) of antibody.
Components:	STAR FLUOR 645 Dye, Cod. F2H085N,DY, Lot NM04A-AR STAR FLUOR Elution buffer, Cod. SEB097,0250, Lot IF27A-NS STAR FLUOR Labeling buffer, Cod. SLB099,0003, Lot IF27B-NS STAR FLUOR Purification Resin, Cod. SPR098, 0250, Lot IF27C-NS Empty plastic column
Storage:	STAR FLUOR Dye at -20°C, other components at +4°C
Stability:	Components provide a kit performance stable for 6 months from date of receipt when stored properly.
Document Release Date:	11/11/2020
Document Revision:	00

Conformance & Quality systems statement

This is to certify that this product conforms to Cyanagen specification.

The product was manufactured via a Quality System certified ISO9001 and tested in accordance with documented quality procedures and approved as a result of meeting the required specification.

Approved:



Quality Control Technician

Reviewed:



Production Manager

Cyanagen warrants that its products conform to the information contained in this and other Cyanagen publications. Purchaser must determine the suitability of the product for its particular use.