

CERTIFICATE OF CONFORMITY

STAR FLUOR 645 Antibody labeling kit

Product Information

Product Number:	F2H085N
Lot Number:	QF13A-XP
Description:	Kit for the fluorescent labelling (Abs 645 nm, Ex 665 nm) of antibody.
Components:	STAR FLUOR 645 Dye, Cod. F2H085N,DY, Lot NM04A-AR STAR FLUOR Elution buffer, Cod. SEB097,0250, Lot OD29B-AR STAR FLUOR Labeling buffer, Cod. SLB099,0003, Lot IF27B-NS STAR FLUOR Purification Resin, Cod. SPR098, 0250, Lot OD29A-AR Empty plastic column
Storage:	STAR FLUOR Dye at -20°C, other components at +4°C
Stability:	Components provide a kit performance stable for 6 months from date of receipt when stored properly.
Document Release Date:	13/06/2023
Document Revision:	00

Conformance & Quality systems statement

This document certifies that the product complies with the manufacturer specifications and its ISO 9001 Quality System.

The product was tested according to documented internal quality procedures and approved since it meets the required specification.

Approved:



Quality Control Technician

Reviewed:



Production Manager

Cyanagen warrants that its products conform to the information contained in this and other Cyanagen publications. Purchaser must determine the suitability of the product for its particular use.

CYANAGEN SRL -

C.F./P.IVA: 02360111203 iscritta nel Registro Imprese di Bologna n. REA: 433637

Capitale sociale sottoscritto e versato € 10.000,00.

"Società soggetta a direzione e coordinamento da parte di Finance & Technology Srl – art.2497 bis cc".