

## CERTIFICATE OF CONFORMITY

### STAR FLUOR 645 Antibody labeling kit

#### Product Information

<b>Product Code:</b>	F2H085N
<b>Lot Number:</b>	QL30B-XP
<b>Description:</b>	Kit for the fluorescent labelling (Abs 645 nm, Ex 665 nm) of antibody.
<b>Components:</b>	STAR FLUOR 645 Dye, Cod. F2H085N,DY, Lot NM04A-AR STAR FLUOR Elution buffer, Cod. SEB097,0250, Lot OD29B-AR STAR FLUOR Labeling buffer, Cod. SLB099,0003, Lot IF27B-NS STAR FLUOR Purification Resin, Cod. SPR098, 0250, Lot OD29A-AR Empty plastic column
<b>Storage:</b>	STAR FLUOR Dye at -20°C, other components at +4°C
<b>Stability:</b>	Components provide a kit performance stable for 6 months from date of receipt when stored properly.
<b>Document Release Date:</b>	30/10/2023
<b>Document Revision:</b>	00

#### Conformance & Quality systems statement

This document certifies that the product complies with the manufacturer specifications and its ISO 9001 Quality System.

The product was tested according to documented internal quality procedures and approved since it meets the required specification.

Approved:   
QC Manager

Reviewed:   
Production Manager

Cyanagen warrants that its products conform to the information contained in this and other Cyanagen publications. Purchaser must determine the suitability of the product for its particular use.  
MO 7.5.1.3 Rev04

