

## CERTIFICATE OF CONFORMITY

### STAR FLUOR MICRO 488 Antibody labeling kit


#### Product Information

<b>Product Code:</b>	F5A122N
<b>Lot Number:</b>	QM09A-XP
<b>Description:</b>	Kit for the fluorescent labelling (Abs 494 nm, Ex 520 nm) of 50-100µg antibody.
<b>Components:</b>	STAR FLUOR MICRO 488 Dye, Cod. F5A122N,DY, Lot QM06A-AR STAR FLUOR 488 Labeling buffer, Cod. SLB147,0003, Lot MB20A-ET STAR FLUOR Purification Spin Column, Cod. SPS121N,001, Lot FJ0793
<b>Storage:</b>	STAR FLUOR MICRO 488 Dye at -20°C, other components at +4°C
<b>Stability/ Expiry Date:</b>	Components provide a kit performance stable up to 6 months from receipt, when properly stored.
<b>Document Release Date:</b>	09/11/2023
<b>Document Revision:</b>	00

#### Conformance & Quality systems statement

This document certifies that the product complies with the manufacturer specifications and its ISO 9001 Quality System.

The product was tested according to documented internal quality procedures and approved since it meets the required specification.

Approved:   
QC Manager

Reviewed:   
Production Manager

Cyanagen warrants that its products conform to the information contained in this and other Cyanagen publications. Purchaser must determine the suitability of the product for its particular use.

MO 7.5.1.3 Rev04