

## CERTIFICATE OF CONFORMITY

### STAR FLUOR MICRO 645 Antibody labeling kit

#### Product Information

<b>Product Number:</b>	F2H124N
<b>Lot Number:</b>	OI06D-XP
<b>Description:</b>	Kit for the fluorescent labelling (Abs 645 nm, Em 665 nm) of 50-100µg antibody.
<b>Components:</b>	STAR FLUOR MICRO 645 Dye, Cod. F2H124N,DY, Lot OH20A-RP STAR FLUOR Labeling buffer, Cod. SLB099,0003, Lot IF27B-NS STAR FLUOR Purification Spin Column, Cod. SPS121N, 001, Lot MF12602 Collection tube
<b>Storage:</b>	STAR FLUOR MICRO 645 – Dye at -20°C, other components at +4°C
<b>Stability/ Expiry Date:</b>	Components provide a kit performance stable up to 6 months from receipt, when properly stored
<b>Document Release Date:</b>	06/09/2021
<b>Document Revision:</b>	00


#### Conformance & Quality systems statement

This document certifies that the product complies with the manufacturer specifications and its ISO 9001 Quality System.

The product was tested according to documented internal quality procedures and approved since it meets the required specification.

Approved: 

Quality Control Technician

Reviewed: 

Production Manager

Cyanagen warrants that its products conform to the information contained in this and other Cyanagen publications. Purchaser must determine the suitability of the product for its particular use.

**CYANAGEN SRL** -

C.F./P.IVA: 02360111203 iscritta nel Registro Imprese di Bologna n. REA: 433637  
Capitale sociale sottoscritto e versato € 10.000,00.

"Società soggetta a direzione e coordinamento da parte di Finance & Technology Srl – art.2497 bis cc".