

## CERTIFICATE OF CONFORMITY

---

### Star Fluor 488 Goat anti-Rabbit

#### Product Information

---

<b>Product Number:</b>	F5A166
<b>Lot Number:</b>	MI10B-FB
<b>Description:</b>	Highly cross adsorbed Goat anti-Rabbit IgG (H&L) - STAR FLUOR 488, conjugated secondary antibody
<b>Components:</b>	Goat anti Rabbit conjugate concentration 1mg/mL in PBS containing 0.01% sodium azide and 10mg/mL BSA after reconstitution
<b>Storage:</b>	+4°C. Protect from light
<b>Stability/ Expiry Date:</b>	10/2022, 6 months after reconstitution at +4°C
<b>Document Release Date:</b>	29/10/2019
<b>Document Revision:</b>	00

#### Conformance & Quality systems statement

---

This is to certify that this product conforms to Cyanagen specification.

The product was manufactured via a Quality System certified ISO9001 and tested in accordance with documented quality procedures and approved as a result of meeting the required specification.

Approved:



Quality Control Technician

Reviewed:



Production Manager

Cyanagen warrants that its products conform to the information contained in this and other Cyanagen publications. Purchaser must determine the suitability of the product for its particular use.