

CERTIFICATE OF CONFORMITY

Star Fluor 645 Goat anti-Rabbit

Product Information

Product Number: F2H168
Lot Number: RP3/227
Description: Highly cross adsorbed Goat anti-Rabbit IgG (H&L) - STAR FLUOR 645, conjugated secondary antibody
Components: Goat anti Rabbit conjugate concentration 1mg/mL in PBS containing 0.01% sodium azide and 10mg/mL BSA after reconstitution
Storage: +4°C. Protect from light
Stability/ Expiry Date: 10/2022, 6 months after reconstitution at +4°C
Document Release Date: 29/10/2019
Document Revision: 00

Conformance & Quality systems statement

This is to certify that this product conforms to Cyanagen specification.

The product was manufactured via a Quality System certified ISO9001 and tested in accordance with documented quality procedures and approved as a result of meeting the required specification.

Approved:


Quality Control Technician

Reviewed:


Production Manager

Cyanagen warrants that its products conform to the information contained in this and other Cyanagen publications. Purchaser must determine the suitability of the product for its particular use.