

## CERTIFICATE OF CONFORMITY

### Star Fluor 680 Goat anti-Rabbit Highly cross-adsorbed

#### Product Information


<b>Product Code:</b>	F2H174
<b>Lot Number:</b>	QM06B-XP
<b>Description:</b>	Highly cross adsorbed Goat anti-Rabbit IgG (H&L) - STAR FLUOR 680, conjugated secondary antibody
<b>Components:</b>	Goat anti Rabbit conjugate concentration 1mg/mL in PBS containing 0.01% sodium azide and 10mg/mL BSA after reconstitution
<b>Storage:</b>	+4°C. Protect from light. Avoid freeze/thaw cycles.
<b>Stability/ Expiry Date:</b>	11/2026, 6 months after reconstitution at +4°C
<b>Document Release Date:</b>	06/11/2023
<b>Document Revision:</b>	00

#### Conformance & Quality systems statement

This document certifies that the product complies with the manufacturer specifications and its ISO 9001 Quality System.

The product was tested according to documented internal quality procedures and approved since it meets the required specification.

Approved:



QC Manager

Reviewed:



Production Manager

Cyanagen warrants that its products conform to the information contained in this and other Cyanagen publications. Purchaser must determine the suitability of the product for its particular use.

MO 7.5.1.3 Rev04