

## CERTIFICATE OF CONFORMITY

### μQPRO-BCA kit Micro

#### Product Information

Product Number:	PRTD2
Lot Number:	PF29A-AR
Description:	Bicinchoninic acid kit for protein determination
Components:	Solution A: Tartrate in alkaline buffer, Lot PF29B-AR Solution B: 4% BCA in water, Lot PF29C-AR Solution C: 4% CuSO <sub>4</sub> 5H <sub>2</sub> O in water, Lot PF29D-AR
Storage:	Room Temperature
Stability/ Expiry Date:	Components provide a kit performance stable up to 12 months from receipt, when properly stored
Document Release Date:	01/07/2022
Document Revision:	00

#### Conformance & Quality systems statement

This document certifies that the product complies with the manufacturer specifications and its ISO 9001 Quality System.

The product was tested according to documented internal quality procedures and approved since it meets the required specification.

Approved:



Quality Control Technician

Reviewed:



Production Manager

Cyanagen warrants that its products conform to the information contained in this and other Cyanagen publications. Purchaser must determine the suitability of the product for its particular use.